

Analytical Performance of the Automated Accelerate Arc™ BC kit and Module for Direct Identification from Positive Blood Cultures using MALDI

Ana Campos-Alvarez¹, Sarah Bolanos¹, Katy Lank¹, Richard Nyberg¹, Shelley, Campeau¹ and Carlos Michel¹

¹Accelerate Diagnostics, Tucson, AZ

Abstract # 01566



BACKGROUND

When treating sepsis, timely identification (ID) is essential to begin targeting appropriate antimicrobial therapy for the causative pathogen. The Accelerate Arc™ BC kit and module is a fully automated solution for cleaning up positive blood culture samples to enable direct identification using matrix-assisted laser desorption ionization (MALDI). This study assessed performance (reportability and accuracy) of the Accelerate Arc™ BC kit (under development) and module for 20 of the most prevalent species isolated from positive blood cultures using MALDI analysis.

METHODS

- A total of 100 isolates (Table 1) were seeded into BD Bactec™ Plus Aerobic blood culture bottles (Becton Dickinson, BD) and incubated on the BD Bactec™ FX blood culture monitoring system.
- After flagging positive, samples were processed on the Accelerate Arc™ BC kit and module (Arc) at 30 min and 60 min time post-positive (TPP) to simulate a typical clinical laboratory workflow.
- 1 µL of Accelerate Arc™-processed sample was spotted in triplicate for MALDI analysis by both the Bruker Biotyper™ (BD) and Vitek® MS (bioMérieux) systems. Gram-positive samples received 1 µL of formic acid prior to addition of matrix, while gram-negative samples received only matrix.
- Growth from the positive blood culture bottle was used to perform conventional MALDI using the Bruker Biotyper™ (e.g., from colony growth) as the comparator.
- Reportability was defined as a sample that had at least 1 of 3 spots that resulted in an ID; if >1 spots gave an ID, they had to agree to be considered reportable (Bruker Biotyper™ = a score of ≥ 1.7; Vitek® MS = any ID).
- Accuracy was defined as ≥1 spot that resulted in an ID for a reportable sample that was in agreement with conventional MALDI ID.

Figure 1. Workflow




Table 1. List of organisms tested

	Species	# Tested
Gram negative	<i>Acinetobacter baumannii</i>	5
	<i>Citrobacter freundii</i>	6
	<i>Citrobacter koseri</i>	4
	<i>Enterobacter cloacae</i> complex	5
	<i>Escherichia coli</i>	6
	<i>Klebsiella aerogenes</i>	5
	<i>Klebsiella oxytoca</i>	5
	<i>Klebsiella pneumoniae</i>	5
	<i>Pseudomonas aeruginosa</i>	6
	<i>Proteus mirabilis</i>	5
	<i>Proteus vulgaris</i>	5
	<i>Serratia marcescens</i>	5
Gram positive	<i>Enterococcus faecalis</i>	5
	<i>Enterococcus faecium</i>	5
	<i>Staphylococcus aureus</i>	6
	<i>Staphylococcus epidermidis</i>	6
	<i>Staphylococcus lugdunensis</i>	4
	<i>Streptococcus agalactiae</i>	4
	<i>Streptococcus pneumoniae</i>	4
	<i>Streptococcus pyogenes</i>	4

RESULTS & CONCLUSIONS

Table 2. Reportability and Accuracy by Organism Group¹

Organism Group	Bruker Biotyper™		Vitek® MS	
	Reportability	Accuracy	Reportability	Accuracy
Gram-negative (n=12 species)	94% (116/124)	100% (116/116)	67% (83/124)	92% (76/83)
Gram-positive (n=8 species)	79% (60/76)	100% (60/60)	71% (54/76)	98% (53/54)
All organisms (n=20 species)	88% (176/200)	100% (176/176)	69% (137/200)	94% (129/137)

¹Data is combined from 30 m & 60 m TPP runs

Table 3. Reportability and Accuracy by Species¹

Organism	Bruker Biotyper™		Vitek® MS	
	Reportability	Accuracy	Reportability	Accuracy
<i>E. coli</i> ²	100%	100%	83%	100%
<i>E. cloacae</i>	100%	100%	90%	89%
<i>K. pneumoniae</i> ²	100%	100%	100%	100%
<i>K. oxytoca</i>	100%	100%	60%	100%
<i>P. mirabilis</i>	100%	100%	90%	100%
<i>P. vulgaris</i>	100%	100%	30%	66%
<i>C. freundii</i>	92%	100%	30%	67%
<i>A. baumannii</i>	90%	100%	40%	100%
<i>K. aerogenes</i>	90%	100%	100%	83%
<i>C. koseri</i>	88%	100%	40%	75%
<i>P. aeruginosa</i> ²	83%	100%	92%	100%
<i>S. marcescens</i>	80%	100%	40%	100%
<i>S. epidermidis</i> ²	100%	100%	100%	100%
<i>S. agalactiae</i>	100%	100%	80%	100%
<i>E. faecium</i> ²	90%	100%	100%	100%
<i>S. aureus</i> ²	83%	100%	83%	100%
<i>E. faecalis</i>	80%	100%	80%	100%
<i>S. pyogenes</i>	75%	100%	13%	100%
<i>S. lugdunensis</i>	50%	100%	38%	100%
<i>S. pneumoniae</i>	35%	100%	50%	67%

¹Data is combined from 30 m & 60 m TPP runs; ²6 most common bloodstream pathogens in original abstract data

- Overall, Arc-processed samples enabled robust reportability and accuracy for 20 bloodstream pathogens on the Bruker Biotyper and Vitek MS systems.
 - Combined reportability across all tested organisms (n=20) was lower for Vitek MS compared to the original data for the 6 most common bloodstream pathogens. This was primarily due to the addition of streptococci and less prevalent gram-negative bacilli (e.g. *P. vulgaris* and *Citrobacter* spp.). Another possibility is instability of some samples post-spotting due to a limitation in the study design where Vitek MS slides were run at an off-site laboratory. However, this would need to be investigated further to rule this out.
 - Out of 200 runs, 2 (1%) had to be repeated due to a technical failure with the instrument.
 - This study was designed to test an equal number of isolates for each of the 20 species evaluated. As a result, overall reportability may vary depending on actual organism prevalence.
- Compared to conventional MALDI, Arc-processed samples had overall agreement of 100% on the Bruker Biotyper and 94% on the Vitek MS.
 - Inaccurate calls on the Vitek MS were either inaccurate to species level but accurate to genus level (e.g. *C. freundii* called as *C. braakii*) or could be mitigated by correlation with gram stain (e.g., *E. cloacae* complex called as *S. capitis*) (data not shown).

ACKNOWLEDGEMENTS

Accelerate Diagnostics, Inc. would like to thank Dr. William Lainhart (Banner University Medical Center, Tucson, AZ) for performing the Vitek MS MALDI runs.

Accelerate Arc™ module is US-IVD and the BC kit is currently US-RUO